






# Navigator

## NTEP Approved Balance



Capacity	610g	6100g
Accuracy Class	III	III
NTEP Number	00-017	00-017
Verification interval (e)	0.1g	1g
Readability (d)	.07g	.7g
Repeatability	0.1g	1g
Linearity (g)	±0.1g	±1g
Weighing modes	g, kg, oz, lb, GN, N, ozt, dwt, ct, tael (3)*, tical*, momme*	
Functions	Parts counting*, percent weighing, checkweighing	
Tare range	to capacity by subtraction	
Stabilization time	3 seconds	
Operating temperature range	50°-104°F / 10°-40°C	
Power requirements	AC Adapter (supplied) or 8- AA batteries (not included)	
Calibration	External calibration mass (not included)	
Display (WxH) (in/mm)	LCD 3.75 x 1.97 / 95 x 50	
Pan size (in/cm)	4.75 / 12 diam.	5.5 x 6 / 14 x 15.2
Dimensions W x H x D (in/cm)	7.25 x 2.25 x 10.5 / 18.4 x 5.7 x 27	
Net Weight (lb/kg)	3.0 / 1.4	3.5 / 1.6
Item No.	<b>N06112</b>	<b>N0F1A2</b>

*\*Available, but not NTEP approved*

*The Navigator models N06112 and N0F1A2 meet or exceed the Class III 6100 division accuracy requirements in accordance with the National Institute of Standards and Technology (NIST) Handbook 44. Certificate of Conformance Number 00-107 was issued under the National Type Evaluation Program (NTEP) of the National Conference of Weights and Measures.*

### Standard Features

RS232 bi-directional interface, auto shut off, weigh below hook, integral security bracket, rugged composite housing, sealed front panel, stainless steel pans, tactile feedback switches, user selectable stability levels and filters, stability indicator, menu and calibration lock-out switch.

### Optional Accessories

Calibration masses, cables, printer, carrying case, security device, scoops. Auxiliary Display (80250787) and In use service cover (80850078).

### Industry Leading Warranty

All Ohaus Navigator balances are protected by a 5 year warranty on labor and parts. Our rugged construction and stringent quality control have been hallmarks of Ohaus scales and balances since 1907.

### Solid Service Support

Should service or repair ever be required, we will perform the necessary work promptly at reasonable rates. In the U.S., use our service hotline 1-800-526-0659.

### Simple To Order

If you would like an Ohaus specialist to assist you in selecting the scale exactly right for you, simply call, in the U.S., 1-800-672-7722... then place your order with the dealer of your choice.

**Ohaus Corporation**  
\* 29 Hanover Road  
Florham Park, NJ 07932 USA  
Tel: 1-800-672-7722  
973-377-9000  
Fax: 973-593-0359

*With offices throughout Europe, Asia and Latin America*

**\* ISO 9001**  
Registered Quality Management System

Internet Address: [www.ohaus.com](http://www.ohaus.com)

LLUS09-00 © Copyright Ohaus Corporation

