



Ruby Scale

Ruby

The Ohaus Ruby Compact Jewelry scale packs powerful features into a convenient portable package. Attractive and versatile with easy operation, higher weighing capacities, long battery life and great performance.

- **Superior RF Protection** – Designed with superior RF protection to deliver accurate and stable results when used in production areas or near cell phones. Approvals for electrical and emissions include CE listed and FCC Class B Part 15 operation.
- **Large Stainless Steel Weighing Platform** – Designed with a large weighing platform to adequately handle larger pieces. In addition, the standard stainless steel platform is removable for easy cleaning.
- **Compact Design** – Low profile design and truly portable when powered by 3 “AA” batteries (included). An optional AC adapter can also be used to power the Ruby scale, when used in a more permanent location.



MODEL	JR300	JR2500
Readability	0.1g	1.0g
Capacity	300g	2500g
Pan Size (in/cm)	5.75 x 5.25/14.6 x 13.3	
Weighing Units	g, lb:oz, dwt, ozt	
Dimensions (in/cm)	8.0 x 5.5 x 1.5 20.3 x 14 x 3.8	

Need pricing or
more information?
OHAUS Customer Service
800.672.7722